TECHNICAL MANUAL


Kartell

H Hydrogenium

N Nitrogenum

C Carbonium

O Oxygenum

P Phosphorum

S Sulphur

Mg Magnesium

Ca Calcium

Na Natrium

K Kalium

Sr Strontium

Ba Barium

Fe Ferrum

Zn Zincum

Cu Cuprum

Pb Plumbum

Ag Argentum

Au Aurum

Pt Platinum

Hg Hydrargirium

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ABBREVIATIONS, TEMPERATURE AND CHEMICAL RESISTANCE OF PLASTICS

The table below lists commonly used abbreviations for plastics. This list covers plastics commonly employed in the manufacture of plastic laboratory ware.

Temperatures appearing in parentheses: limits tolerated for intervals only.

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LIST OF PLASTICS AND THEIR CHEMICAL RESISTANCES TO SUBSTANCE GROUPS

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“CLEAN” AND ECOLOGICAL PLASTICS

Thanks to technological innovation and progressive adaptation to the most elevated eco-compatibility standards, Kartell plastic raw material are of excellent quality. In fact the used plastic materials grant wide compatibilities; hereunder you can find some standards conformities related to their characteristics.

- Foodstuff contact suitable (National and European Directive CE 1935/2004)
- Foodstuff contact suitable (US FDA CFR 170/199)
- PTFE: foodstuff contact suitable (FDA CFR TITLE 177.1550)
- Absence or limitation according to the International Directives of: heavy metals, biphenyls and phthalates
- RoHS Directives Conformity
- BSE and/or TSE Directives Conformity
- Atex Directives Conformity

Kartell Labware promotes the research in the field of the plastic materials, through the direct contact with the most important worldwide raw materials manufacturers and the research of innovative materials. We would like to underline that Kartell Labware materials, if not chemically contaminated, are totally recyclable

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STERILISING PLASTIC LABORATORY WARE

Before sterilising any items of plastic laboratory ware, verify that no contamination or residues are present. Their presence could destroy plastics during sterilisation or autoclaving.

Observe the tolerated temperature ranges of plastic when autoclaving plastic laboratory ware. Remove any stoppers, fittings, or caps from plastic laboratory ware prior to autoclaving. Plastic vessels should be autoclaved separately from their closures and other fitting. Autoclaving vessels with their closures in place can lead to deformation and destruction of the vessels. All statements are advisory only, and imply no liability on the part of Kartell.

All statements relating to the resistances of plastic laboratory ware to high temperatures, chemicals, and to sterilisation and cleaning procedures have been cautiously formulated, based on statements of raw material manufactures, on statements appearing in the literature, plus experience gained in actual practice.

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CLEANING PLASTIC LABORATORY WARE

All polyolefins, such as LDPE, HDPE, PP and PMP (TPX), as well as the fluorinated hydrocarbons PTFE, PFA, FEP, ETFE and E-CTFE, have wettable surfaces that are both highly resistant to high temperatures and chemical attack and easy to clean. Slight contamination can be removed using a chemically neutral (pH 7) cleaning agent. Heavy contamination can be removed using an alkaline (pH up to 12) cleaning agent.

Use only chemically neutral (pH7) cleaning agents on polycarbonate (PC) or Polystyrene (PS).

Laboratory ware used in trace analyses should be cleaned in a 1-N hydrochloric acid (HCl) solution for periods of not more than 6 hours and then rinsed in distilled water in order to preclude contamination by cations or anions.

Never use scouring powders or abrasive sponges when cleaning plastic laboratory ware. Use no alkaline cleaning agents on polycarbonate (PC) laboratory ware.

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